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
## Introduction

The Stability Trainer is used to improve balance, posture, and coordination. Training on the unstable surface requires balance and a strong posture. It is mainly used in prevention (children and seniors), in rehabilitation after an injury, and in sports in order to enhance performance.

Studies prove that targeted balance training can reduce the frequency of falls and secondary injuries. Due to higher demands placed on supporting muscles, the balance trainer can also improve your performance with conventional exercises.

## Stability Trainer

### Materials and characteristics

The colour-coded foam pads include an anti-slip coating and are available in three levels of firmness. The material is kind to the skin (PVC-free), light-weight and enables an important base for a healthy training. The footprint of 40 x 50 cm allows each kind of individual coordination- and stabilisation exercises.

Material	Skill level	
firm	easy	
medium	medium	
soft	hard	

### Precautions

Every time you train, first verify that the Stability Trainer is resting firmly on the surface below it and cannot slide. Avoid training on slippery floors or rugs. The Stability Trainer may be placed onto a gymnastics mat in order to prevent slipping. However, this does increase the level of difficulty for exercises. Persons suffering from pronounced balance, posture, or sensory disturbances and persons with high risk of falls (caused by medications, osteoporosis, or reduced muscle strength) should only train under the instruction of a physician, physiotherapist, or specialised trainer.

- » Remove sharp or pointed objects from the immediate vicinity of the training area.
- » Regular cleaning with warm water and soap is recommended. For disinfectant use the ARTZT disinfectant spray. The Stability Trainer can be allowed to air-dry. Avoid using rough or high-friction cloths or sponges when cleaning the Stability Trainer.
- » The Stability Trainer should be stored flat and without any items on top of it.
- » Avoid exposure to direct sunlight when storing the Stability Trainer.
- » You should train barefoot or with socks, as shoes might damage the Stability Trainer.

### Training tips

For effective training to improve and enhance your balance, progress from basic to advanced exercises:

- » Use a sturdy, fixed object for support while performing the exercises (for example, you can hold on to a wall or a sturdy chair).
- » Train close to the chosen support object so that you can hold yourself up if necessary.
- » Progress to train without any support at all.
- » Progress from a firm surface to the Stability Trainers.

Exercises with the Stability Trainer may also be completed at increasing levels of difficulty:

- Level 1: Firm surface without Stability Trainer
- Level 2: Green Stability Trainer
- Level 3: Blue Stability Trainer
- Level 4: Black Stability Trainer

Intermediate levels can be adapted by stacking Stability Trainers from one colour.

The increase in training intensity should be tailored to individual performance so that the exercises are always performed safely, correctly, and with good posture and coordination.

### Training tips

Health-oriented, preventive training focuses on the improvement of the general sense of balance. Beginners should select 6 to 8 of the following exercises in order to achieve this goal. Static balance exercises (exercises without movement) should be held 3 times for 10-30 seconds at a time. Dynamic exercises (leg swinging exercises, stepping, jumping, etc.) should be repeated 10-15 times in 1-3 sets. You should never experience EXTREMELY heavy strain while training. Once you have established a foundation and the exercises come easily, you can increase the level of difficulty, number of repetitions (sets), and training frequency. Please observe the following principles:

- » Maintain correct posture and conduct movements in a slow, controlled manner in the initial position, during the exercise, and at the end.
- » Training should be pain-free. Take a break if an exercise becomes painful. You should consult a healthcare provider if pain occurs repeatedly.
- » Take short breaks or change sides between static exercises. You may also alternate between dynamic and static exercises.
- » Examine each exercise in detail and visualise the process.
- » ALWAYS exercise both the right and left side of the body even when the illustration only shows one side.

**Übung - Untere Extremität | Exercises for the lower extremities**  
**Parallelstand | Stand in parallel position**

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Find a stable standing position (without wobbling). Look straight ahead.

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Keep your position stable and look down at your feet. Change your line of vision: Straight ahead, at your feet, straight ahead...

**Parallelstand | Stand in parallel position**

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Shift your weight from the left to the right leg and vice versa.

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Keep your position stable and close your eyes. Now shift your weight in different directions while keeping your eyes closed.

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Stand on tiptoe and then bend your knees slightly. Return to the initial position.

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Stand on a Stability Trainer with both feet. Bend down so that your fingertips touch the floor.

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Stand on tiptoe from this bent over position. The arms are raised over the head. Your eyes should follow your hands.

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Put your front foot onto the Stability Trainer (slightly bend the knee). Shift your weight back and forth between the front and back legs.

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Briefly hold the front and back positions, respectively.

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Shift your weight forward and slightly raise the heel of the free leg.

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Shift your weight forward and lift the foot of the free leg off the floor.

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Shift your weight to the supporting leg and move the free leg forward and backward slightly.

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Shift your weight to the supporting leg and move the free leg back and forth.

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Shift your weight to the supporting leg, move the other leg forward and briefly touch the heel to the floor (your weight remains on the supporting leg!). Bring your front leg back.

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Exercises a-g may also be performed with the arm insert. For sports training, exercises a-g may also be performed on tiptoe (supporting leg).

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Place the back foot onto the Stability Trainer. Keep the supporting leg firmly on the floor. Push off with the toes of the back foot.

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Stand approximately 50 cm in front of the Stability Trainer.

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Now take a large flat-footed step onto the Stability Trainer and back.

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ATTENTION: Be sure to use a non-slip surface so that the Stability Trainer does not slide during this dynamic exercise.

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Jump onto or over the Stability Trainer. Once again, be sure to use a non-slip surface!

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While standing on one leg, move the right hand to the knee of the opposite leg.

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Perform opposing arm and leg movements while standing on one leg.

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Place a Stability Trainer under both hands. You can use a towel or a mat for cushioning your knees.

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Find a stable position and perform push-ups.

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Place a Stability Trainer under both hands. Perform the following movements:

- » Shift your weight from the left to the right.
- » Alternately lift your hands.
- » Perform push-ups.

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Get down on all fours with a Stability Trainer under one knee. Perform the following movements:

- » Extend and flex the free leg.
- » Extend and flex the free leg and opposite arm.
- » Same as above, but with an additional Stability Trainer underneath the supporting hand.

# UPGRADE YOUR VITALITY!

## Produktprogramm ARTZT vitality® | Product range ARTZT vitality® Medical Fitness - Made in Germany

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH.  
Unauthorized use is strictly prohibited.  
© 2009 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

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- Flowin® Personal Training
- Gymstick™ Health
- BOSU® Balance Training

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ARTZT vitality® ist ein eingetragenes Warenzeichen der Ludwig Artzt GmbH.
© 2009 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.
der Ludwig Artzt GmbH.